



Recombinant Human T-cell leukemia virus 2 Protein Rex

Product Code	CSB-YP768043HQL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q85601
Product Type	Recombinant Protein
Immunogen Species	Human T-cell leukemia virus 2 (HTLV-2)
Purity	>85% (SDS-PAGE)
Sequence	MPKTRRQRTR RARRNRPPTP WPISQDLDR SYMDTPSTCL AIVYRPIGVP SQVVYVPPAY IDMPSWPPVQ STNSPGTSM DAL SALLSNT LSLASPPSP REPQGSRSL PLPPLLSPPR FHLPSFNQCE STPPTMDAW NQPSGISSPP SPSPNLASVP KTSTPPGKPK
Source	Yeast
Protein Names	Recommended name: Protein Rex Alternative name(s): Rev homolog Rex-2
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.