



Recombinant Human Basic immunoglobulin-like variable motif-containing protein (BIVM)

Product Code	CSB-MP768208HU
Abbreviation	BIVM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UB2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPNVAETERS NDSNGGEHKS ERKSPEENLQ GAVKSFCTSA SGAPLGPKGD GHYPWSCPVT HTREKIYAIC SDYAFLNQAT SIYKTPNPSR SPCLPDSTSL SAGNNSSRYI GIPTSTSEII YNEENLENL SNSLGKLPLA WEIDKSEFDG VTTNSKHKSG NAKKQVSKRK TSDKKGRYQK ECPQHSPLED IKQRKVLDLR RWYICISRPQY KTSCGISSLI SCWNFLYSTM GAGNLPPITQ EEALHILGFQ PPFEDIRFGP FTGNTTLMRW FRQINDHFHV KGCSYVLYKP HGKNKTAGET ASGALSKLTR GLKDESLAYI YHCQNHYFCP IGFEATPVKA NKAFSRGP PQEVEYWILI GESSRKHPAI HCKKWADIVT DLNTQNPEYL DIRHLERGLQ YRKTCKVGGN LHCIIAFQRL NWQRFGLWNF PFGTIRQESQ PPTHAQGIK SESEDNISKK QHGRLGRSFS ASFHQDSAWK KMSSIHERRN SGYQGYSDYD GND
Source	Mammalian cell
Target Names	BIVM
Protein Names	Recommended name: Basic immunoglobulin-like variable motif-containing protein
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.