



Recombinant Human Serpin A9 (SERPINA9)

Product Code	CSB-YP768223HU
Abbreviation	SERPINA9
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86WD7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ANAPSAY PRPSSTKSTP ASQVYSLNTD FAFRLYRRLV LETPSQNIFF SPVSVSTSLA MLSLGAHSVT KTQILQGLGF NLTHTPESAI HQGFQHLVHS LTVPSKDLTL KMGSAFVKK ELQLQANFLG NVKRLYEAEV FSTDFSNPSI AQARINSHVK KKTQGKVVDI IQGLDLLTAM VLVNHIFKKA KWEKPFHPEY TRKNFPFLVG EQVTVHVPMM HQKEQFAFGV DTELNCFVLQ MDYKGDVAVF FVLPSKGKMR QLEQALSART LRKWSHSLQK RWIEVFIPRF SISASYNLET ILPKMGIQNV FDKNADFSGI AKRDSLQVSK ATHKAVLDVS EEGTEATAAT TTKFIVRSKD GPSYFTVSFN RTFLMMITNK ATDGILFLGK VENPTKS
Source	Yeast
Target Names	SERPINA9
Protein Names	Recommended name: Serpin A9 Alternative name(s): Centerin Germinal center B-cell-expressed transcript 1 protein
Expression Region	24-417
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.