



Recombinant Human Uncharacterized protein CXorf67 (CXorf67)

Product Code	CSB-MP768227HU
Abbreviation	CXorf67
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86X51
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MATQSDMEKE QKHQQDEGQG GLNNETALAS GDACGTGNQD PAASVTTVSS QASPSGGAAL SSSTAGSSAA AATSAAIFIT DEASGLPIIA AVLTERHSDR QDCRSPHEVF GCVVPEGGSQ AAVGPQKATG HADEHLAQTK SPGNSRRRKQ PCRNQAAPAQ KPPGRRLFPE PLPPSSPGFR PSSYPCSGAS TSSQATQPGP ALLSHASEAR PATRSRITLV ASALRRRASG PGPVIRGCTA QPGPAFPHRA THLDPARLSP ESAPGPARRG RASVPGPARR GCDSAPGPAR RGRDSAPVSA PRGRDSAPGS ARRGRDSAPG PALRVRTARS DAGHRSTSTT PGTGLRSRST QQRSALLSRR SLSGSADENP SCGTGSERLA FQSRSGSPDP EVPSRASPPV WHAVRMRASS PSPPGRFFLP IPQQWDESSS SSYASNSSSP SRSPGLSPSS PSPEFLGLRS ISTPSPESLR YALMPEFYAL SPVPPEEQAE IESTAHPATP PEP
Source	Mammalian cell
Target Names	CXorf67
Protein Names	Recommended name: Uncharacterized protein CXorf67
Expression Region	1-503
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.