



Recombinant Human Nitrilase homolog 1 (NIT1)

Product Code	CSB-MP768228HU
Abbreviation	NIT1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q86X76
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLGFITRPPH RFLSLLCPGL RIPQLSVLCA QPRPRAMAISSSSCELPLVA VCQVTSTPDK QQNFKTCAEL VREAARLGAC LAFLPEAFDF IARDPAETLH LSEPLGGKLL EEYTQLAREC GLWLSLGGFH ERGQDWEQTQ KIYNCHVLLN SKGAVVATYR KTHLCDVEIP GQGPMCESNS TMPGPSLESP VSTPAGKIGL AVCYDMRFPE LSLALAQAGA EILTYPSTAFG SITGPAHWEV LLRARAIETQ CYVAAAQCG RHHEKRASYG HSMVVDPWGT VVARCSEGGP LCLARIDLNY LRQLRRHLPV FQHRRPDLYG NLGHPLS
Source	Mammalian cell
Target Names	NIT1
Protein Names	Recommended name: Nitrilase homolog 1 EC= 3.5.-.-
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>