



Recombinant Human AN1-type zinc finger protein 4 (ZFAND4)

Product Code	CSB-EP768229HU-B
Abbreviation	ZFAND4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86XD8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDNRKEPPFF NDDNMGPFYY RLHFCDTMEL FIETLTGTCTF ELRVSPFETV ISVKAKIRRL EGIPICRQHL IWNNMELEND YCLNDYNISE GCTLKLVLAM RGGPINTRRV PTDDPLRKMA EYLDSSRVEV WEKTSCSKQV TFLVYQEGDQ LNFFPAVDRG DGLTPLSDS SKKIDFHLHV LRRKGEHRMS GGSMYNSDTD EDEETEPSSS GQQIENSIT MNKMKLLKAK MKNMNLKSKP KKAVKIKPHP PVAPRPSSGS TAPSRHLLR VLPNIGQSCS PAFGNAYPPE ISRNGISSLA TQLSAERYIS SITGEFLKED NSWENNTLSH FSSNVKLPQ IPHLELGNDQ ELADSVLHLG SSLPRQTKHF LGNLPSSNGN IVLPSEECVT EQSLLPKVGS LASFAEGNAD EQSSGLEGAC KVNLELLTN ADKGLKAPEQ HLKHAVAGVLN GESVETSVLN YRELSPHKNR LLSPLRCSAP MSLHNSLVKP ERQSKCFEFG KLQPSSSQSL DVQNITDSSF SRTTCFQGVK VDSLGRSDV ISKVEARDIT EMTNKASKEP VGCVNNISFL ASLAGSTSRN RLQSTRGAGR LQNSGTGLST NLQHFQEENF RKSSPQLEHT GVFLSTHGVG MNGNNAAGK SVGECTTHHL PPVKAPLQTK KKTTHCFLC GKKTGLASSY ECRCGNNFCA SHRYAETHGC TYDYKSAGR YLHEANPVVN APKLPKI
Source	E.coli
Target Names	ZFAND4
Protein Names	Recommended name: AN1-type zinc finger protein 4 Alternative name(s): AN1-type zinc finger and ubiquitin domain-containing protein-like 1
Expression Region	1-727
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.