



# Recombinant Human AN1-type zinc finger protein 4 (ZFAND4)

<b>Product Code</b>	CSB-MP768229HU
<b>Abbreviation</b>	ZFAND4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86XD8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MDNRKEPPFF NDDNMGPFYY RLHFCDTMEL FIETLTGTCTF ELRVSPFETV ISVKAKIRRL EGIPICRQHL IWNNMELEND YCLNDYNISE GCTLKLVLAM RGGPINTRRV PTDDPLRKMA EYLDSSRVEV WEKTSCSKQV TFLVYQEGDQ LNFFPAVDRG DGLTPLSDS SKKIDFHLHV LRRKGEHRMS GGSMYNSDTD EDEETEPSSS GQQIENSIT MNKMLLKAK MKNMNLKSKP KKAVKIKPHP PVAPRPSSGS TAPSRHLLR VLPNIGQSCS PAFGNAYPPE ISRNGISSLA TQLSAERYIS SITGEFLKED NSWENNTLSH FSSNVKLPQ IPHLELGNDQ ELADSVLHLG SSLPRQTKHF LGNLPSSNGN IVLPSEECVT EQSLLPKVGS LASFAEGNAD EQSSGLEGAC KVNLELLTN ADKGLKAPEQ HLKHVAGVLN GESVETSVLN YRELSPHKNR LLSPLRCSAP MSLHNSLVKP ERQSKCFEFG KLQPSSSQSL DVQNITDSSF SRTTCFQGVK VDSLGRSDV ISKVEARDIT EMTNKASKEP VGCVNNISFL ASLAGSTRN RLQSTRGAGR LQNSGTGLST NLQHFQEENF RKSSPQLEHT GVFLSTHGVG MNGNNAAGK SVGECTTHHL PPVKAPLQTK KKTTHCFLC GKKTGLASSY ECRCGNNFCA SHRYAETHGC TYDYKSAGRR YLHEANPVVN APKLPKI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZFAND4
<b>Protein Names</b>	Recommended name: AN1-type zinc finger protein 4 Alternative name(s): AN1-type zinc finger and ubiquitin domain-containing protein-like 1
<b>Expression Region</b>	1-727
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.