



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-B2-1 (CYCB2-1)

Product Code	CSB-EP768696OFG-B
Abbreviation	CYCB2-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XSJ6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDRASENRRRL AAVGKPVPGI GEMGNRRPLR DINNLVGAPP HPSAIKKPM LEKSGKEEQK PALVVSHRPM TRNFAASLTR KEQLDHQVSV ADAAVVCTDP QKNPIPDGTV DDDVESCESN DYIAVDECND TDEDESMMDI DSADSGNPLA ATEYVEELYK FYRENEEMSC VQPDYMSSQG DINEKMRAIL IDWLIEVHHK FELMDETLFL TVNIVDRFLE KQVVPKRLQ LVGVTAMLLA CKYEEVAVPV VEDLVLISDR AYTGGQILEM EKLILNTLQF NMSVPTPYVF MRRFLKAAQS DKQLQLLSFF ILELSLVEYQ MLKYRPSLLA AAVYTAQCA LTRCQQWTKT CELHSRYTGE QLLECSRMMV DFHQKAGAGK LTGVHRKYST FKFGCAAQTE PALFLLESGA GGYNLQKQPC
Source	E.coli
Target Names	CYCB2-1
Protein Names	Recommended name: Cyclin-B2-1 Alternative name(s): CycB2-Os1 G2/mitotic-specific cyclin-B2-1 Short name= CycB2;1
Expression Region	1-420
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.