



Recombinant Human Ankyrin repeat domain-containing protein 33 (ANKRD33)

Product Code	CSB-MP768761HU
Abbreviation	ANKRD33
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z3H0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVACYHGFQS VVALLSHCPF LDVNQQDKGG DTALMLAAQA GHVPLVLLLL NYYVGLDLER RDQRGLTALM KAAMRNRCEC VATLLMAGAD LTAVDPVRGK TALEWAVLTD SFDTVWRIRQ LLRRPQVEQL SQHYKPEWPA LSGLVAAQAQA QAQVAPSLLE RLQATLSLPF APSPQEGGVL DHLVTATTSL ASPFVTTACH TLCPDHPPSL GTRSKSVPEL LVPAAEQSFR TPKSGPSSLA IPGAQDREEE TGGGGQNGTE VGEDGIGQAG NR
Source	Mammalian cell
Target Names	ANKRD33
Protein Names	Recommended name: Ankyrin repeat domain-containing protein 33
Expression Region	1-272
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.