



Recombinant Human Neuroepithelial cell-transforming gene 1 protein (NET1)

Product Code	CSB-BP768775HU
Abbreviation	NET1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z628
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEPELAAQKQ PRPRRRSRRA SGLSTEGATG PSADTSGSEL DGRCSLRRGS SFTFLTPGPN WDFTLKRKRR EKDDDVVLSL SLDLKEPSNK RVRPLARVTS LANLISPVRN GAVRRFGQTI QSFTLRGDHR SPASAQKFSS RSTVPTPAKR RSSALWSEML DITMKESLTT REIRRQEAIY EMSRGEQDLI EDLKLARKAY HDPMLKLSIM SEEELTHIFG DLDSYIPLHE DLLTRIGEAT KPDGTVEQIG HILVSWLPRL NAYRGYCSNQ LAAKALLDQK KQDPRVQDFL QRCLESPFSR KLDLWSFLDI PRSRLVKYPL LLKEILKHTP KEHPDVQLLE DAILIIQGVL SDINLKKGES ECQYYIDKLE YLDEKQRDPR IEASKVLLCH GELRSKSGHK LYIFLFQDIL VLTRPVTRNE RHSYQVYRQP IPVQELVLED LQDGDVRMGG SFRGAFSNSE KAKNIFRIRF HDPSPAQSHT LQANDVFHKQ QWFNCIRAAI APFQSAGSPP ELQGLPELHE ECEGNHPSAR KLTAQRRAST VSSVTQVEVD ENAYRCGSGM QMAEDSKSLK THQTQPGIRR ARDKALSGGK RKETLV
Source	Baculovirus
Target Names	NET1
Protein Names	Recommended name: Neuroepithelial cell-transforming gene 1 protein Alternative name(s): Proto-oncogene p65 Net1 Rho guanine nucleotide exchange factor 8
Expression Region	1-596
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.