



# Recombinant Human Neuroepithelial cell-transforming gene 1 protein (NET1)

<b>Product Code</b>	CSB-YP768775HU
<b>Abbreviation</b>	NET1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z628
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEPELAAQKQ PRPRRSRRA SGLSTEGATG PSADTSGSEL DGRCSLRGGS SFTFLTPGPN WDFTLKRKR EKDDVVVLS SLDLKEPSNK RVRPLARVTS LANLISPVRN GAVRRFGQTI QSFTLRGDHR SPASAQKFSS RSTVPTPAKR RSSALWSEML DITMKESLTT REIRRQEAIY EMSRGEQDLI EDLKLARKAY HDPMLKLSIM SEEELTHIFG DLDSYIPLHE DLLTRIGEAT KPDGTVEQIG HILVSWLPRL NAYRGYCSNQ LAAKALLDQK KQDPRVQDFL QRCLESPFSR KLDLWSFLDI PRSRLVKYPL LLKEILKHTP KEHPDVQLLE DAILIIQGVL SDINLKKGES ECQYYIDKLE YLDEKQRDPR IEASKVLLCH GELRSKSGHK LYIFLFQDIL VLTRPVTRNE RHSYQVYRQP IPVQELVLED LQDGDVRMGG SFRGAFSNSE KAKNIFRIRF HDPSPAQSHT LQANDVFHKQ QWFNCIRAAI APFQSAGSPP ELQGLPELHE ECEGNHPSAR KLTAQRRAST VSSVTQVEVD ENAYRCGSGM QMAEDSKSLK THQTQPGIRR ARDKALSGGK RKETLV
<b>Source</b>	Yeast
<b>Target Names</b>	NET1
<b>Protein Names</b>	Recommended name: Neuroepithelial cell-transforming gene 1 protein Alternative name(s): Proto-oncogene p65 Net1 Rho guanine nucleotide exchange factor 8
<b>Expression Region</b>	1-596
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.