



# Recombinant Danio rerio Uncharacterized protein C7orf26 homolog (zgc:56409)

<b>Product Code</b>	CSB-MP768795DIL
<b>Abbreviation</b>	zgc:56409
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7ZTY6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDIRHALLR RDPLSAAKEV LYHLDISLGS ALQSSAGPPA GLEKSTVELV EEFIFHVSKD RNIQSKRMSC VQELQLLEII CSYFQEQSKD AIRQVIFSAL FSLQGNKADE SRMAMLGKLV SMSIAVCRVP ILECAAIWLQ RTHSAWCVRL AQVLMEDYCT LVPCAISTLQ NICSASPRFC CQLITAVTAL YDFSTDELTP PPALLEMLVG WITEDPRLLL LTFINSPLNS SLPLGCLEIT PLLGLLRWCV KAPLAYHRGR KSAMTNGHGD GEKGTAYEEL YSKLHLSVLQ VFLMLQVHLT EKNMFCRLTV LPVESVAALI EEVGRLCCKL APLLDASHIQ LALDRLAQUAL QVAMATGALL CAREDLRALC SRLPHNLLQ LVMSGPVVQP PPHSSPFQPN MYPHIHPGRP SPLSPHSPHQ STLSSPHSPH TVLTAHPHP ALAPHRPLTA HATHPALSPH AFHPASIAFP YRPIR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	zgc:56409
<b>Protein Names</b>	Recommended name: Uncharacterized protein C7orf26 homolog
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.