



Recombinant *Xenopus laevis* Target of rapamycin complex 2 subunit MAPKAP1 (mapkap1), partial

Product Code	CSB-MP768827XBE
Abbreviation	mapkap1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZYH1
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Source	Mammalian cell
Target Names	mapkap1
Protein Names	Recommended name: Target of rapamycin complex 2 subunit MAPKAP1 Short name= TORC2 subunit MAPKAP1 Alternative name(s): Mitogen-activated protein kinase 2-associated protein 1 Stress-activated map kinase-interacting protein 1 Short n
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.