



Recombinant Danio rerio DDB1- and CUL4-associated factor 13 (dcaf13)

Product Code	CSB-MP768843DIL
Abbreviation	dcaf13
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q803X4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKVKVLCRNP DDYVRETRD IQRVPRNYDP TLHPFEVPRE YTRALNATKL ERVFAKPFVA SLDGHRDVIS CITKHAKSLS TVISGACDGE VKVWNLPKRE CTR TVQAHEG FVRGICSRFC GTSFFTVGDD KTIKQWNMEA PGYGVREEPI NTILGKAVFT GIDHHQREGT FVTCGQTVDI WDEQRSSPVL SFSWGVDSFS SVRYNPVEPD LLASCASDRS IVLYDTRESA PLRKVIMQLR SNTLCWNPME AYYFTCANED YNLYTYDIRN LDVPVTVHMD HVSALVDVDY SPTGREFVSA SFDKTIRIFP KDKGHSREVV HTKRMQHVIC VKWSADSKFI MSGSDEMNI LWKANASEKL GVLSTREKTA ANYNKKLIQK FQHHPQVRR ARHRHLPRDV LKQKRELKEM KEARRRKEQN VRKHSKPGSV PLLTEKEKHV VKVVE
Source	Mammalian cell
Target Names	dcaf13
Protein Names	Recommended name: DDB1- and CUL4-associated factor 13 Alternative name(s): WD repeat and SOF domain-containing protein 1
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.