



Recombinant Mouse Signal transduction protein CBL-C (Cblc)

Product Code	CSB-MP768908MO
Abbreviation	Cblc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q80XL1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MAAAAAPRGW QRGEPRALSR AVKLLQRLEE QCRDPRMVTG PPSLRDLLPR TAQLLGEVAK ARREAREDPE GPGGADDFLA IYLANLEVKG RQVAELLPPR GKKDQVNDVF REGSRFRRQL AKLALIFSHM HAELSALFPA GKCYGHLQYL TKGSAHIFWR QNCGVRCVLP WAEFQSLCA CHPVEPGPTM QALRSTDLT CNGHVSVEFE DVFTRLFQPW PTLLRNWQLL AVNHPGYMAF LTYDEVQTRL QAYRDKPGSY IFRPSCTRLG QWAIGYVSSD GSILQTIPLN KPLLQVLLKG QKDGIFLFPD GKKHNPDLTE LCRVEPYQRI QVSEEQLLLY QAMNSTFQLC KICAERDKDV RIEPCGHLLC SCCLAAWQDS DSQTCPCFCRC EIKGREAVSI CQAQERPTEV RTAADGSRDN CHQEAAEQKL GPVIPSAPSL LPEDQFPQGP QDKGWLT LAP LALPRLRPPL PLPKMASVLW EVTSRPRARE EATESS
Source	Mammalian cell
Target Names	Cblc
Protein Names	Recommended name: Signal transduction protein CBL-C Alternative name(s): SH3-binding protein CBL-3 SH3-binding protein CBL-C
Expression Region	1-496
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.