




Recombinant Mouse Tenascin (Tnc), partial

Product Code	CSB-YP768917MO
Relevance	Extracellular matrix protein implicated in guidance of migrating neurons as well as axons during development, synaptic plasticity as well as neuronal regeneration. Promotes neurite outgrowth when provided to neurons in culture. May play a role in supporting the growth of epithelial tumors. Ligand for integrins ITGA8:ITGB1, ITGA9:ITGB1, ITGAV:ITGB3 and ITGAV:ITGB6. In tumors, stimulates angiogenesis by elongation, migration and sprouting of endothelial cells
Abbreviation	Recombinant Mouse Tnc protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q80YX1
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	GLLYPFPRDCSQAMLNGDTTSGLYTIYINGDKTQALEVYCDMTSDGGGWIVFL RRKNGREDFYRNWKAYAAGFGDRREEFWLGLDNLSKITAQQQYELRVLDLQD HGESAYAVYDRFSVGDASKSRYKLVKVEGYSGTAGDSMNYHNGRSFSTYDKDT DSAITNCALSYKGAFWYKNCHRNVNLMGRYGDNNHSQGVNWFHWKGHEYSIQ FAEMKLRPSN
Research Area	Signal Transduction
Source	Yeast
Target Names	Tnc
Protein Names	Hexabrachion (Tenascin-C) (TN-C) (Hxb)
Expression Region	1884-2099aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	26.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.