



Recombinant UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase (murE)

Product Code	CSB-BP769471SZB
Abbreviation	murE
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83MG0
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MADRNLRDLL APWVPDAPSR ALREMTLDSR VAAAGDLFVA VVGHQADGRR YIPQAIAGGV AAIIEAKDE ATDGEIREMH GVPVIYLSQL NERLSALAGR FYHEPSDNLR LVGVTGTNGK TTTTQLLAQW SQLLGETSAV MGTVGNLLG KVIPTENTTG SAVDVQHELA GLVDQGATFC AMEVSSHGLV QHRVAALKFP ASVFTNLSRD HLDYHGDMEH YEAAKWLLYS EHHCGQAIIN ADDEVGRRWL AKLPDAVAVS MEDHINPNCH GRWLKATEVN YHDSGATIRF SSSWGDGEIE SHLMGAFNVS NLLLALATLL ALGYPLADLL KTAARLQPVC GRMEVFTAPG KPTVVVDYAH TPDALAKALQ AARLHCAGKL WCVFGCGGDR DKGKRPLMGA IAEEFADVAV VTDDNPRTEE PRAIINDILA GMLDAGHAKV MEGRAEAVTC AVMQAKENDV VLVAGKGHED YQIVGNQRLD YSDRVTVARL LGGIA
Source	Baculovirus
Target Names	murE
Protein Names	Recommended name: UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase EC= 6.3.2.13 Alternative name(s): Meso-A2pm-adding enzyme Meso-diaminopimelate-adding enzyme UDP-MurNAc-L-Ala-D-Glu:meso-diaminopimelate ligase
Expression Region	1-495
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.