



Recombinant Arabidopsis thaliana Nuclear transcription factor Y subunit A-7 (NFYA7)

Product Code	CSB-MP769551DOA
Abbreviation	NFYA7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84JP1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MTSSIHESLD NIGSHEKQEQ RDSHFQPPIP SARNYESIVT SLVYSDPGTT NSMAPGQYPY PDPYYRSIFA PPPQPYTGVB LQLMGVQQQG VPLPSDAVEE PVFVNAKQYH GILRRRQSRA RLESQNKVIK SRKPYLHESR HLHAIRPRG CGGRFLNAKK EDEHHEDSSH EEKSNLSAGK SAMAASSGTS
Source	Mammalian cell
Target Names	NFYA7
Protein Names	Recommended name: Nuclear transcription factor Y subunit A-7 Short name= AtNF-YA-7
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.