



Recombinant *Oryza sativa* subsp. japonica PHD finger protein ALFIN-LIKE 3 (Os03g0818300, LOC_Os03g60390)

Product Code	CSB-MP769574OFG
Abbreviation	Os03g0818300, LOC_Os03g60390
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84TV4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEMAPAAQVA SNPRTVEDIF KDYSARRGAL VRALTSDVDE FFGLCDPDKE NLCLYGLANG SWEVALPAEE VPPELPEPAL GINFARDGMN RRDWLSLVAV HSDSWLVSA FFFAARLNGN ERKRLFNMIN DLPTVYEALV DRKHVRDRSG VDSSGKSKHS TKRTGEGQVK RSRVVAEEYE DDDEEHNETF CGTCGGLYNA NEFWIGCDIC ERWFHGKCVR ITPAKAEHIK HYKCPDCSSS SSKKTRL
Source	Mammalian cell
Target Names	Os03g0818300
Protein Names	Recommended name: PHD finger protein ALFIN-LIKE 3
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.