



# Recombinant Human Pre-mRNA-processing factor 39 (PRPF39)

<b>Product Code</b>	CSB-EP769774HU
<b>Abbreviation</b>	PRPF39
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86UA1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQNSHMDEYR NSSNGSTGNS SEVVVEHPTD FSTEIMNVTE MEQSPDDSPN VNASTEETEM ASAVDLPVTL TETEANFPPE YEKFWKTVEN NPQDFTGWVY LLQYVEQENH LMAARKAFDR FFIHYPYCYG YWKKYADLEK RHDNIKPSDE VYRRGLQAIP LSVDLWIHYI NFLKETLDPG DPETNNTIRG TFEHAVLAAG TDFRSDRLWE MYINWENEQG NLREVTAIYD RILGIPTQLY SHHFQRFKEH VQNNLPRDLL TGEQFIQLRR ELASVNGHSG DDGPPGDDL SGIEDITDPA KLITEIENMR HRIIEHQEM FNYNEHEVSK RWTFEEGIKR PYFHVKPLEK AQLKNWKEYL EFEIENGTHE RVVVLFERCV ISCALYEEFW IKYAKYMENH SIEGVRHVFS RACTIHLPKK PMVHMLWAAF EEQQGNINEA RNILKTFEEC VLGLAMVRLR RVSLERRHGN LEEAEHLLQD AIKNAKSNE SSFYAVKLAR HLFKIQKNLP KSRKVLLEAI ERDKENTKLY LNLLEMEYS DLKQNEENIL NCFDKAVHGS LPIKMRITFS QRKVEFLEDF GSDVNKLLNA YDEHQTLLE QDSLKRKAEN GSEPEEKKA HTEDTTSSST QMIDGDLQAN QAVYNYSAWY QYNYQNPWNY GQYYPPPPT
<b>Source</b>	E.coli
<b>Target Names</b>	PRPF39
<b>Protein Names</b>	Recommended name: Pre-mRNA-processing factor 39 Alternative name(s): PRP39 homolog
<b>Expression Region</b>	1-669
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.