



Recombinant Human Pre-mRNA-processing factor 39 (PRPF39)

Product Code	CSB-MP769774HU
Abbreviation	PRPF39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UA1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQNSHMDEYR NSSNGSTGNS SEVVVEHPTD FSTEIMNVTE MEQSPDDSPN VNASTEETEM ASAVDLPVTL TETEANFPPE YEKFWKTVEN NPQDFTGWVY LLQYVEQENH LMAARKAFDR FFIHYPYCYG YWKKYADLEK RHDNIKPSDE VYRRGLQAIP LSVDLWIHYI NFLKETLDPG DPETNNTIRG TFEHAVLAAG TDFRSDRLWE MYINWENEQG NLREVTAIYD RILGIPTQLY SHHFQRFKEH VQNNLPRDLL TGEQFIQLRR ELASVNGHSG DDGPPGDDL SGIEDITDPA KLITEIENMR HRIIEIHQEM FNYNEHEVSK RWTFEEGIKR PYFHVKPLEK AQLKNWKEYL EFEIENGTHE RVVVLFERCV ISCALYEEFW IKYAKYMENH SIEGVRHVFS RACTIHLPKK PMVHMLWAAF EEQQGNINEA RNILKTFEEC VLGLAMVRLR RVSLERRHGN LEEAEHLLQD AIKNAKSNE SSFYAVKLAR HLFKIQKNLP KSRKVLLEAI ERDKENTKLY LNLLEMEYSG DLKQNEENIL NCFDKAVHGS LPIKMRITFS QRKVEFLEDF GSDVNKLLNA YDEHQTLLEKE QDSLKRKAEN GSEEPKKA HTEDTTSSST QMIDGDLQAN QAVYNYSAWY QYNYQNPWNY GQYYPPPPT</p>
Source	Mammalian cell
Target Names	PRPF39
Protein Names	Recommended name: Pre-mRNA-processing factor 39 Alternative name(s): PRP39 homolog
Expression Region	1-669
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.