



Recombinant Human Radial spoke head protein 3 homolog (RSPH3)

Product Code	CSB-EP769775HU-B
Abbreviation	RSPH3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86UC2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTVKPAKAAS LARNLAKRRR TYLGGAAGRS QEPEVPCAAV LPGKPGDRNC PEFPPDRTL GCWATDAAPA AGLCGAGSEP SIAPTSCAGN LPSRPPPLLS PLLASRNPCP WHYHLHLSGSH NTLAPTCKFA KLHRKRGSQP PDMASALTDR TSRAPSTYTY TSRPRALPCQ RSRYRDSLQ PDEEPMHYGN IMYDRRVIRG NTYALQTGPL LGRPDSLELQ RQREARKRAL ARKQAQEQLR PQTPEPVEGR KHVDVQTELY LEEIADRIIE VDMECQTD AF LDRPPTPLFI PAKTGKDVAT QILEGELDFD DLEV KPVLEV LVGKTIEQSL LEVMEEEELA NLRASQREYE ELRNSERA EV QRLEE QERRH REEKERRKKQ QWEIMHKHNE TSQKIAARAF AQRYLADLLP SVFGSLRDSG YFYDPIERDI EIGFLPWLMN EVEKTMEYSM VGRTVLDMLI REVVEKRLCM YEHGEDTHQS PEPEDEPGGP GAMTESLEAS EFLEQSMSQT RELLLDGGYL QRTTYDRRSS QERKFMEERE LLGQDEETAM RKSLGEEELS
Source	E.coli
Target Names	RSPH3
Protein Names	Recommended name: Radial spoke head protein 3 homolog Alternative name(s): Radial spoke head-like protein 2
Expression Region	1-560
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.