



# Recombinant Human Radial spoke head protein 3 homolog (RSPH3)

<b>Product Code</b>	CSB-MP769775HU
<b>Abbreviation</b>	RSPH3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q86UC2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTVKPAKAAS LARNLAKRRR TYLGGAAGRS QEPEVPCAAV LPGKPGDRNC PEFPPDRTL GCWATDAAPA AGLCGAGSEP SIAPTSCAGN LPSRPPPLLS PLLASRNPCP WHYHLHLSGSH NTLAPTCFKA KLHRKRGSQP PDMASALTDR TSRAPSTYTY TSRPRALPCQ RSRYRDSLQ PDEEPMHYGN IMYDRRVIRG NTYALQTGPL LGRPDSLELQ RQREARKRAL ARKQAQEQLR PQTPEPVEGR KHVDVQTELY LEEIADRIIE VDMECQTDAF LDRPPTPLFI PAKTGKDVAT QILEGELDFD DLEVKPVLEV LVGKTIEQSL LEVMEEEELA NLRASQREYE ELRNSERA EV QRLEEQERRH REEKERRKKQ QWEIMHKHNE TSQKIAARAF AQRYLADLLP SVFGSLRDSG YFYDPIERDI EIGFLPWLMN EVEKTMEYSM VGRTVLDMLI REVVEKRLCM YEHGEDTHQS PEPEDEPGGP GAMTESLEAS EFLEQSMSQT RELLLDGGYL QRTTYDRRSS QERKFMEERE LLGQDEETAM RKSLGEEELS
<b>Source</b>	Mammalian cell
<b>Target Names</b>	RSPH3
<b>Protein Names</b>	Recommended name: Radial spoke head protein 3 homolog Alternative name(s): Radial spoke head-like protein 2
<b>Expression Region</b>	1-560
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.