



Recombinant Human Exocyst complex component 3-like protein (EXOC3L1)

Product Code	CSB-MP769787HU
Abbreviation	EXOC3L1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VI1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSAAKDEMQ PALSPGPEWP EQERAEQLAR GAALKWASGI FYRPEQLARL GQYRSREVQR TCSLESRLKS VMQSYLEGVQ TGVWQLAQAI EVVQGTREAL SQARGLLQGM SQALQTLEPL RERVAQHKQL QALSHLLPRL RAVPAAVSHT QTLIDGQQFL EAYVSLRELE QLREDTWAPL GGLELPVFQG LDLLFEALGQ AVEAAAGAAG KLAREDPALL VAAVRVAEVE TGRTTPLGQV PRDWRQRCLR ALQEGLEQAH FGSPLLPAPG ALPGWLEALR VALPVELATA EALVAPCCPP QYNVVQLWAH TLHSGLRRLSL QNLLAGPELE AADAFALLHW ALHVYLGQEM MGSLELGPEA DVSQLEPLLT LENIEQLEAT FVANIQASVS QWLQNALDGE VAEWGREHGP NTDPSPGSYYS PMPAIVLQIL EENIRVASLV SESLQQRVHG MALSELGTFL RSFSDALIRF SRDHFRGKSM APHYVPYLLA ALNHKSALSS SVSVLQLDGA PSGALAPVEA ALDELQRRYI RLVLEALQAE LQPLFADLPS RQWLSSPELL QSVCERTGRF CRDFWRVRNP TVQLLLAEAE RAVVLQYLSA LMQGRLVCRG ADERTQAAER LRHDAAQLQQ LFLSLGLEEN AHCAPVLLAL RELLNLRDPA LLGLEVAGLR QQFPDVSEDH VSALLGLRGD LSREQHLAAL SSLQAALPPS PRASRRVLFS LVPAPALAPA SCLPSGSCAR ALLLAE
Source	Mammalian cell
Target Names	EXOC3L1
Protein Names	Recommended name: Exocyst complex component 3-like protein Alternative name(s): Protein Jiangli
Expression Region	1-746
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.