



Recombinant Human Lebercilin (LCA5)

Product Code	CSB-BP769789HU
Abbreviation	LCA5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VQ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p> MGERAGSPGT DQERKAGKHH YSYLSDFETP QSSGRSSLVS SSPASVRRKN PKRQTSDGQV HHQAPRKPS KGLPNRKGVR VGFRSQSLNR EPLRKDIDL TKRILSARLL KINELQNEVS ELQVKLAELL KENKSLKRLQ YRQEKALNKF EDAENEISQL IFRHNNEITA LKERLRKSQE KERATEKRVK DTESELFRTK FSLQKLKEIS EARHLPERDD LAKKLVS AEL KLDDTERRIK ELSKNLELST NSFQRQLLAE RKRAYEAHDE NKVLQKEVQR LYHKLKEKER ELDIKNIYSN RLPKSSPNKE KELALRKNA CQSDFADLCT KGVQTMEDFK PEEYPLTPET IMCYENKWEE PGHLTLDLQS QKQDRHGEAG ILNPIMEREE KFVTDEELHV VKQEVKLED EWEREELDKK QKEKASLLER EEKPEWETGR YQLGMYP IQN MDKLQGE EEE RLKREMLLAK LNEIDRELQD SRNLKYPVLP LLPDFESKLH SPERSPKTYR FSESSERLFN GHHLQDISFS TPKGEGQNSG NVRSPASPNE FAFGSYVPSF AKTSERSNPF SQKSSFLDFQ RNSMEKLSKD GVDLITRKEK KANLMEQLFG ASGSSTISSK SSDPNSVASS KGDIDPLNFL PGNKGSRDQE HDEDEGFFLS EGRSFNPNRH RLKHADDKPA VKAADSVEDE IEEVALR </p>
Source	Baculovirus
Target Names	LCA5
Protein Names	Recommended name: Lebercilin Alternative name(s): Leber congenital amaurosis 5 protein
Expression Region	1-697
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This gene encodes a protein that is thought to be involved in centrosomal or ciliary functions. Mutations in this gene cause Leber congenital amaurosis type V. Alternatively spliced transcript variants are described.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.