



Recombinant Human Protein Hook homolog 3 (HOOK3)

Product Code	CSB-BP769790HU
Abbreviation	HOOK3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VS8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MFSVESLERA ELCESLLTWI QTFNVDAPCQ TVEDLTNGVV MAQVLQKIDP AYFDENWLNR IKTEVGDNWR LKISNLKKIL KGILDYNHEI LGQQINDFTL PDVNLIGEHS DAAELGRMLQ LILGCAVNCE QKQEYIQAIM MMEESVQHVV MTAIQELMSK ESPVSAGNDA YVDLDRQLKK TTEELNEALS AKEEIAQRCH ELDMQVAALQ EEKSSLLAEN QVLMERLNQS DSIEDPNSPA GRRHLQLQTQ LEQLQEETFR LEAAKDDYRI RCEELEKEIS ELRQQNDEL TLADEAQLSK DEIDVLRHSS DKVSKLEGQV ESYKKKLEDL GDLRRQVKLL EEKNTMYMQN TVSLEELRK ANAARSQLET YKRQVVELQN RLSEESKAD KLD FEYKRLK EKVDSLQKEK DRLRTERDSL KETIEELRCV QAQEGQLTTQ GLMPLGSQES SDSLAAEIVT PEIREKLIRL QHENKMLKLN QEGSDNEKIA LLQSLLDDAN LRKNELETEN RLVNQRLLLEV QSQVEELQKS LQDQGSKAED SVLLKKKLEE HLEKLHEANN ELQKKRAIIE DLEPRFNSS LKIEELQEAL RKKEEEMKQM EERYKKYLEK AKSVIRTLDP KQNQGAPEI QALKNQLQER DRLFHSLEKE YEKTKSQREM EEKYIVSAWY NMGMTLHKKA AEDRLASTGS GQSFLARQRQ ATSSRRSYPG HVQPATAR
Source	Baculovirus
Target Names	HOOK3
Protein Names	Recommended name: Protein Hook homolog 3 Short name= h-hook3 Short name= hHK3
Expression Region	1-718
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.