



Recombinant Human THO complex subunit 6 homolog (THOC6)

Product Code	CSB-EP769793HU
Abbreviation	THOC6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86W42
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MERAVPLAVP LGQTEVFQAL QRLHMTIFSQ SVSPCGKFLA AGNNYGQIAI FSLSSALSSE AKEESKPPVV TFQAHDGPVY SMVSTDRHLL SAGDGEVKAW LWAEMLKKGK KELWRRQPPY RTSLEVPEIN ALLLVPKENS LILAGGDCQL HTMDLETGTF TRVLRGHTDY IHCLALRERS PEVLGGEDG AVRLWDLRTA KEVQTIEVYK HEECSRPHNG RWIGCLATDS DWMVCGGGPA LTLWHLRSST PTTIFPIRAP QKHVTFYQDL ILSAGQGRCV NQWQLSGELK AQVPGSSPGL LSLSLNQQPA APECKVLTAA GNSCRVDVFT NLGYRAFSL S F
Source	E.coli
Target Names	THOC6
Protein Names	Recommended name: THO complex subunit 6 homolog Alternative name(s): Functional spliceosome-associated protein 35 Short name= fSAP35 WD repeat-containing protein 58
Expression Region	1-341
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.