



Recombinant Human Prostate tumor-overexpressed gene 1 protein (PTOV1)

Product Code	CSB-MP769809HU
Abbreviation	PTOV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YD1
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MVRPRRAPYR SGAGGPLGGR GRPPRPLVVR AVRSRSWPAS PRGPQPPRIR ARSAPPMEGA RVFGALGPIG PSSPGLTLGG LAVSEHRLSN KLLAWSGVLE WQEKRRPYSD STAKLKRTLQ CQAYVNQGEN LETDQWPQKL IMQLIPQQLL TTLGPLFRNS QLAQFHFTNR DCDSLKGLCR IMGNGFAGCM LFPHISPCEV RVLMLLYSSK KKIFMGLIPY DQSGFVSAIR QVITTRKQAV GPGGVNSGPV QIVNNKFLAW SGVMEWQEPR PEPNSRSKRW LPSHVYVNQG EILRTEQWPR KLYMQLIPQQ LLTTLVPLFR NSRLVQFHFT KDLETLSLC RIMDNFGAGC VHFSYKASCE IRVLMMLLYSS EKKIFIGLIP HDQGNFVNGI RRVIANQQQV LQRNLEQEQQ QRGMGG
Source	Mammalian cell
Target Names	PTOV1
Protein Names	Recommended name: Prostate tumor-overexpressed gene 1 protein Short name= PTOV-1 Alternative name(s): Activator interaction domain-containing protein 2
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.