



Recombinant Human Myosin light chain kinase family member 4 (MYLK4)

Product Code	CSB-EP769814HU-B
Abbreviation	MYLK4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86YV6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLKVKRLEEF NTCYNSNQLK KMAFFQCREE VEKVKCFLEK NSGDQDSRSG HNEAKEVWSN ADLTERMPVK SKRTSALAVD IPAPPAPFDH RIVTAKQGAV NSFYTVSKTE ILGGGRFGQV HKCEETATGL KLAAKIIKTR GMKDKEEVKN EISVMNQLDH ANLIQLYDAF ESKNDIVLVM EYVDGGELFD RIIDESYNLT ELDTILFMKQ ICEGIRHMHQ MYILHLDLKP ENILCVNRDA KQIKIDFGL ARRYPREKL KVNFGTPEFL APEVVNYDFV SFPTDMWSVG VIAYMLLSGL SPFLGDNDAAE TLNNILACRW DLEDEEFQDI SEEAKEFISK LLIKEKSWRI SASEALKHPW LSDHKLHSRL NAQKKKNRGS DAQDFVTK
Source	E.coli
Target Names	MYLK4
Protein Names	Recommended name: Myosin light chain kinase family member 4 EC= 2.7.11.1 Alternative name(s): Sugen kinase 85 Short name= SgK085
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.