



# Recombinant Mouse X-ray repair cross-complementing protein 5 (Xrcc5)

<b>Product Code</b>	CSB-EP026233MO-B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P27641
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAWSGNKAAVVLCDVGVAMGNSFPGEEESPIEQAKKVMTMFVQRQVFSESK DEIALVLYGTDGTDNALAGKDQYQNITVCRHMLPDLLEDIGNKIQPSSQQA DFLDALIVCMDLIQRETIGKKFGKKHIEVFTDLSSPFSQDQLDVIICNLKKGISL QFFLPFPIDKNGEPGERGDLDSGLDHLKPSFPQKGLTEQQKEGIRMVTRVMLS LEGEDGLDEIYSFSESLRQLCVFKKIERRSMPWPCQLTIGPNLSIKIVAYKSIVQ EKFKKSWVVVDARTLKKEDIQKETVYCLNDDDETEVSKEDTIQGYRYGSDIIPF SKVDEEQMKYKSEGKCFVGLGFCSSQVHRRFFMGHQVLKVFAAKDDEAAA VALSSLVHALDELNMVAIVRYAYDKRSNPQVGVAFPYIKDAYECLVYVQLPFM EDLRQYMFSSLKNNKCTPTEAQLSAIDDLIDMSLVKKNEEEDIVEDLFPTSKI PNPEFQRLYQCLLHRALHLQERLPPIQQHILNMLDPPTMKAACESPLSKVKTL FPLTEVIKKKNQVTAQDVFQDNHEEGPAAKKYKTEKEEDHISISLAEGNITKV GSVNPVENFRFLVRQKIASFEEASLQLISHIEQFLDTNETLYFMKSMDCIKAFRE EAIQFSEEQRFSNFLEALREKVEIKQLNHFWIIVVQDGVTLITKDEGPGSSITAE EATKFLAPKDKAKEDTTGPEEAGDVDDLLDMI
<b>Source</b>	E.coli
<b>Target Names</b>	Xrcc5
<b>Protein Names</b>	Recommended name: X-ray repair cross-complementing protein 5 EC= 3.6.4.- Alternative name(s): ATP-dependent DNA helicase 2 subunit 2 ATP-dependent DNA helicase II 80 kDa subunit CTC box-binding factor 85 kDa subunit Short name=
<b>Expression Region</b>	1-732aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.