



Recombinant *Oryza sativa* subsp. japonica Expansin-B17 (EXPB17)

Product Code	CSB-BP770259OFG
Abbreviation	EXPB17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7X6J9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	SFLF DGGKSKSAAA AAAMDMEWRP ATATWYGDAE GDGSTGGACG YGSLVDVPM KARVGSVSPV LFKDGEGCGA CYKVKCLDHG ICSRRAVTVI VTDECPGGLC AFGRTHFDLS GAAFSRMAVA GAGGHLRDRG QLSVVYRRTA CKYGGKNI AF RVNEGSTNFW LLLVEFEDG QGDIGSMQIK QANSVEWLDM KHVWGATWCL VRGPLVGPFS VRLTTLSAQK ALTARDVIPR NWKPTATYTS RLNFEAAL
Source	Baculovirus
Target Names	EXPB17
Protein Names	Recommended name: Expansin-B17 Alternative name(s): Beta-expansin-17 OsEXPB17 OsaEXPb1.13
Expression Region	27-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.