



Recombinant *Oryza sativa* subsp. japonica Clathrin light chain 1 (Os04g0679100, LOC_Os04g58240)

Product Code	CSB-YP770274OFG
Abbreviation	Os04g0679100, LOC_Os04g58240
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XKE9
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASFFADDGA DELPRTASHP FDADDDAAPD ASGGAAADD T GYGGYASFVD GGVEDVEEEE EEIAVESEGV PIGHVSGGFS PSPFSPDPEL DGGDGPILPP PAQMGAEEGI LLREWRRQNA IVLEEKERKE KELRAQILAE AEEFKKAFYE KRIQNCETNK VHNREKIF VAGQEKFAE ADKQYWKSIS ELIPHEIATI EKRGKKDKDK KPSITVIQGP KPGKPTDLR MRQILVKLKH APPPHMMQPP PASAADGAK DGAADGTPAP ANGTKKPAES KEKPANGSPA EAEKEQPAAS E
Source	Yeast
Target Names	Os04g0679100
Protein Names	Recommended name: Clathrin light chain 1
Expression Region	1-301
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.