



# Recombinant *Oryza sativa* subsp. japonica Probable gamma-aminobutyrate transaminase 2, mitochondrial (Os04g0614500, LOC\_Os04g52440)

<b>Product Code</b>	CSB-YP770277OFG
<b>Abbreviation</b>	Os04g0614500, LOC_Os04g52440
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7XN12
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. japonica (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	DST EGNQFKGHSM LAPFTAGWHS TDLEPLIIEE SEGSYVYDSK GNKYLDLTLG LWCTALGGSE PRLVKAATDQ LNKLPFYHSF WNSTAKPPLD LAEELISMFT AKEMGKVFFT NSGSEANDSQ VKLVWYYNNA LGRPNKKKII AQSQAYHGST LISASLSGLP AMHLKFDLPA PFVLHTDCPH YWRFGLPGEA EEEFATRLAD NLENLILKEG PETVAAFAIE PVIGAGGVIP PPKTYFEKIQ AVLQKYDVLV IADEVITGFG RLGTMFGSDL YNIKPDVSL AKALSSAYVP IGATLVSPEI SDVVHSQSNK IGFFAHGFTY SGHPVSCAVA LEALKIYRER NIPAHVKQIS PRFQEGIKAF AGSSIIGETR GVGLLLATEF ANNKSPNDPF PVEWGVAQIF GAECKKRGML VKVVGDEIAM SPPLIMSQRE VDGLVSIYGE ALKATEERVA ELRSKKK
<b>Source</b>	Yeast
<b>Target Names</b>	Os04g0614500
<b>Protein Names</b>	Recommended name: Probable gamma-aminobutyrate transaminase 2, mitochondrial EC= 2.6.1.-
<b>Expression Region</b>	38-497
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.