



Recombinant Human HEAT repeat-containing protein 3 (HEATR3)

Product Code	CSB-MP770350HU
Abbreviation	HEATR3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z4Q2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGKSRTKRFK RPQFSPTGDC QAEAAAAANG TGGEEDDGPA AELLEKLQHP SAEVRECACA GLARLVQQRP ALPGLARRDA VRRLGPLLLD PSLAVRETAA GALRNLSACG GFEVCDDMVT KDIMTPLVAL LKECSAGLDS NEMSLQEKKD QNRNSIENIA NETVNVLWNI CECSSRAVSI FNKEGCLEIV LKYLRSRFTN VDLAISVAYC LQTVTEDNPE LLKSFSATAL NMLESALLSP VSSMESLLLK TLVAGTIWNL KDIIPCKSQA EIINALLKIL SEVLGMDAGE MVIQMKEAET QRLKTAAEAE EILENTNGDD LIEDEMEGI SHKRRVRRKT FVSDLLPPTD KELRETIALL TAQQTALEII VNMCCNEDPS DDEWEELSSS DESDAFMENS FSECGQLFS PLCLSHEVHT ALTNYLIPKK IFEKTAFPNS IAVDLCSRNP TWKPLIRKMN TIQCRAFCL QSLVSLLDVE HLGGAALQT LAQHLSQLLF SQPDFAKHVD FLEAISSALR ALLQTMASKN ISQCMTDPQL MTLCKAGIHS SNVGVRVNVV SILGITGSVL AKEDGTLET LKNIGCFLEEV TTKDPSLVVA GEALDALFDV FADGKEAERA SIQIKLLSAL KEFQPVFKMK IRKEGRGNYS TDQLCVLDNV KMNLRRFIAY QETVEKRLTS
Source	Mammalian cell
Target Names	HEATR3
Protein Names	Recommended name: HEAT repeat-containing protein 3
Expression Region	1-680
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.