



# Recombinant Human HEAT repeat-containing protein 3 (HEATR3)

<b>Product Code</b>	CSB-YP770350HU
<b>Abbreviation</b>	HEATR3
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z4Q2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGKSRTKRFK RPQFSPTGDC QAEAAAAANG TGGEEDDGPA AELLEKLQHP SAEVRECACA GLARLVQQRP ALPGLARRDA VRRLGPLLLD PSLAVRETAA GALRNLSACG GFEVCDDMVT KDIMTPLVAL LKECSAGLDS NEMSLQEKKD QNRNSIENIA NETVNVLWNI CECSSRAVSI FNKEGCLEIV LKYLRSRPTN VDLAISVAYC LQTVTEDNPE LLKSFSATAL NMLESALLSP VSSMESLLLK TLVAGTIWNL KDIIPCKSQA EIINALLKIL SEVLGMDAGE MVIQMKEAET QRLKTAEEAE EILENTNGDD LIEDEMEGI SHKRRVRRKT FVSDLLPPTD KELRETIALL TAQQTALIEI VNMCCNEDPS DDEWEELSSS DESDAFMENS FSECGQLFS PLCLSHEVHT ALTNYLIPKK IFEKTAFPNS IAVDLCSRNP TWKPLIRKMN TIQCRAFCL QSLVSLLDVE HLGGAALQT LAQHLSQLLF SQPDFAKHVD FLEAISSALR ALLQTMASKN ISQCMTDPQL MTLCKAGIHS SNVGVRVNVV SILGITGSVL AKEDGTLET LKNIGCFLEEV TTKDPSLVVA GEALDALFDV FADGKEAERA SIQIKLLSAL KEFQPVFKMK IRKEGRGNYS TDQLCVLDNV KMNLRRFIAY QETVEKRLTS
<b>Source</b>	Yeast
<b>Target Names</b>	HEATR3
<b>Protein Names</b>	Recommended name: HEAT repeat-containing protein 3
<b>Expression Region</b>	1-680
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.