



# Recombinant Human Hydroxysteroid 11-beta-dehydrogenase 1-like protein (HSD11B1L)

<b>Product Code</b>	CSB-EP770354HU-B
<b>Abbreviation</b>	HSD11B1L
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z5J1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	YYWDD NFDPASLQGA RVLLTGANAG VGEELAYHYA RLGSHLVLTA HTEALLQKVV GNCRKLGAPK VFYIAADMAS PEAPESVVQF ALDKLGGLDY LVLNHIGGAP AGTRARSPQA TRWLMQVNFV SYVQLTSRAL PSLTDSKGS LVVSSLLGRV PTSFSTPYSA AKFALDGGFFG SLRRELDVQD VNVAITMCVL GLRDRASAAE AVRSTSRPR QPEHRGVPLQ SQTAMFLPPT VPGARTLTET PLRGWPQPKM KSSRQKSKTE KNDGHLEPVT AWEVQVPRVR RLCRGLARPH LFGHD
<b>Source</b>	E.coli
<b>Target Names</b>	HSD11B1L
<b>Protein Names</b>	Recommended name: Hydroxysteroid 11-beta-dehydrogenase 1-like protein EC= 1.1.1.- Alternative name(s): 11-beta-hydroxysteroid dehydrogenase type 3 Short name= 11-DH3 Short name= 11-beta-HSD3 Short-chain dehydrogenase/reduc
<b>Expression Region</b>	16-315
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.