



Recombinant Human Spermatogenesis-associated protein 19, mitochondrial (SPATA19)

Product Code	CSB-MP770355HU
Abbreviation	SPATA19
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z5L4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SSDIDV VESEAVSVLH HWLKKTEEEA SRGIKEKLSI NHPSQGVREK MSTDSPPTHG QDIHVTRDVV KHHLSKSDLL ANQSQEVLEE RTRIQFIRWS HTRIFQVPSE MTEDIMRDRI EQVRRSISRL TDVSAQDFSM RPSSSDC
Source	Mammalian cell
Target Names	SPATA19
Protein Names	Recommended name: Spermatogenesis-associated protein 19, mitochondrial Alternative name(s): Spermatogenic cell-specific gene 1 protein Short name= Spergen-1
Expression Region	25-167
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.