



Recombinant Human Sorting nexin-20 (SNX20)

Product Code	CSB-EP770357HU-B
Abbreviation	SNX20
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q7Z614
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MASPEHPGSP GCMGPITQCT ARTQQEAPAT GPDLPHPGPD GHLDTHSGLS SNSSMTTREL QQYWQNQKCR WKHVKLLFEI ASARIEERKV SKFVVYQIIV IQTGSFDNNK AVLERRYSD F AKLQKALLKT FREEIEDVEF PRKHLTG NFA EEMICERRRA LQEYLGLLYA IRCVRRSREF LDFLTRPELR EAFGCLRAGQ YPRALELLLR VLPLQEKLTA HCPAAAVPAL CAVLLCHRD L DRPAEFAAAG ERALQRLQAR EGHRYYAPLL DAMVRLAYAL GKDFVTLQER LEESQLRRPT PRGITLKELT VREYLH
Source	E.coli
Target Names	SNX20
Protein Names	Recommended name: Sorting nexin-20 Alternative name(s): Selectin ligand-interactor cytoplasmic 1
Expression Region	1-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>