



Recombinant Human UPF0469 protein KIAA0907 (KIAA0907)

Product Code	CSB-EP770364HU-B
Abbreviation	KIAA0907
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z7F0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSAGSATHPG AGGRRSKWDQ PAPAPLLFLP PAAPGGEVTS SGGSPGGTTA APSGALDAAA AVAAKINAML MAKGKPKPTQ NASEKLQAPG KGLTSNKS KD DLVVAEVEIN DVPLTCRNL TRGQTQDEIS RLSGAAVSTR GRFMTTEEKA KVGPGDRPLY LHVQGQTREL VDRAVNRIKE IITNGVVKAA TGTSPTFNGA TVTVYHQPAP IAQLSPAVSQ KPPFQSGMHY VQDKLFGVLE HAVPTFNVKE KVEGPGCSYL QHIQIETGAK VFLRGKGGSGC IEPASGREAF EPMYIYISHP KPEGLAAAKK LCENLLQTVH AEYSRFVNQI NTAVPLPGYT QPSAISSVPP QPPYYPSNGY QSGYPVPPP QQPVPQPPYGV PSIVPPAVSL APGVLPALPT GVPPVPTQYP ITQVQPPAST GQSPMGGPFI PAAPVKTALP AGPQPQPQPQ PPLPSQPQAQ KRRFTEELPD ERESGLLG YQ HGPIHMTNLG TGFSSQNEIE GAGSKPASSS GKERERDRQL MPPPAFPVTG IKTESDERNG SGTLTGSHDY PAKKMKTTEK GFGLVAYAAD SSDEEEHGG HKNASSFPQG WSLG YQYPSS QPRAKQQMPF WMAP
Source	E.coli
Target Names	KHDC4
Protein Names	Recommended name: UPF0469 protein KIAA0907
Expression Region	1-614
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.