



Recombinant *Saccharomyces mikatae* Dihydroorotate dehydrogenase (fumarate) (URA1)

Product Code	CSB-MP770367SBAK
Abbreviation	URA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z893
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces mikatae</i> (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTASLTTKFL NNTYENPFMN ASGVHCMTTE ELDELADSKA GAFITKSATT LEREGNPKPR YISVPLGSIN SMGLPNEGID YYLSYVLNRQ KKYPDAPAIF FSVAGMSIDE NLNLLKKIQD SEFNGITELN LSCPNVPGKP QVAYDFELTK ETLEKVFVFF KKPLGIKLPP YFDFAHFDII AKILNEFPLA YVNSINSIGN GLFIDVEKES VVVKPKNGFG GIGGEYVKPT ALANVRAFYT RLRPDIKVIQ TGGIKSGKDA FEHLLCGASM LQIGTELQKE GVKIFERIEK ELKDIMEAKG YTSIDQFRGK LNSL
Source	Mammalian cell
Target Names	URA1
Protein Names	Recommended name: Dihydroorotate dehydrogenase (fumarate) Short name= DHOD Short name= DHODase Short name= DHODEHASE EC= 1.3.98.1 Alternative name(s): Dihydroorotate oxidase
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.