



Recombinant Danio rerio Pancreas transcription factor 1 subunit alpha (ptf1a)

| | |
|--------------------------|---|
| Product Code | CSB-EP770375DIL |
| Abbreviation | ptf1a |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q7ZSX3 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MDTVLDPFTG LDSFSSSYFD DDDFFTDHSS RDHLDTDDFL DEDVDLFTNQ IQEYYKDSRI SQDGDYCDVG NFSFSSSSST FSYGCADSTS ELSPHRDGG LKRRRRMRSE VEMQQLRQAA NVRERRRMQS INDAFEGLRS HIPTLPYEKR LSKVDTLRLA IGYINFLAEL VQSDMPIRNP HSDALNQPKK VIICHRGTRS PSPNDPDYGL PPLAGHSLSW TDEKQLKDQN IIRTAKVWTP EDPRKLHLKS SINNIENEPP FNFIS |
| Source | E.coli |
| Target Names | ptf1a |
| Protein Names | Recommended name: Pancreas transcription factor 1 subunit alpha Alternative name(s): Pancreas-specific transcription factor 1a bHLH transcription factor p48 |
| Expression Region | 1-265 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |