



Recombinant Danio rerio Inositol-tetrakisphosphate 1-kinase (itpk1)

Product Code	CSB-MP770379DIL
Abbreviation	itpk1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZU91
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MQTFVKGKRV GYWLSEKKIK KLNFTQFVDL CRKQGIEMIQ LDLSQPIESQ GPFDVIIHKL TDHIVDADQN VTESLLLQVQ VQDYIDAHPE TVILDPLPAI RTLLDRCKSY KLIHLEHSM EDDRICSPPF MVLKTECGFE TLEQLHKHGI TFPFICKPQV AHGTNSHEMA IIFSEEDLKD IKPPCVLQSF INHNAVLYKV FVVGAEYSVV QRPSIRNFPS GPTDRRAISF NSHHVSKPES SSHLTCRDNM VGQSWKPSNE VIQKISRKLH QALGISLFGI DIIINNQTGQ HAVIDINAFP GYEGVPEFFD DLLSHISSVL QGQVCNGVAC GHLRVNGTAQ SQTVHCGMLG NDSSWLMDSE GIKKGP HQGL SCCGGC MTPN FHQHTRSSLT AETSSQ
Source	Mammalian cell
Target Names	itpk1
Protein Names	Recommended name: Inositol-tetrakisphosphate 1-kinase EC= 2.7.1.134 Alternative name(s): Inositol 1,3,4-trisphosphate 5/6-kinase Short name= Inositol-trisphosphate 5/6-kinase Short name= Ins(1,3,4)P(3) 5/6-kinase EC= 2.7
Expression Region	1-396
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.