



Recombinant Danio rerio E3 ubiquitin-protein ligase RNF181 (rnf181)

Product Code	CSB-MP770393DIL
Abbreviation	rnf181
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZW78
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MASYFDEHDC EPTNPEEQYR QNALLELARS LMQGLDIDSG SFDLSDWDQR LPPPAAKAVV QSLPVVIISP EQADKGVKCP VCLLEFEEQE SVREMPCKHL FHTGCILPWL NKTNSCPLCR LELPTDNADY EEFKKDKERR RQREHRLEDL HGAMYT
Source	Mammalian cell
Target Names	rnf181
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF181 EC= 6.3.2.- Alternative name(s): RING finger protein 181
Expression Region	1-156
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.