



Recombinant Danio rerio Nuclear distribution protein nudE-like 1-B (ndel1b)

Product Code	CSB-MP770423DIL
Abbreviation	ndel1b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q803Q2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MDTEMMPKFS SKDEEIDYWK CLSLKYKKSC HDAQEELQEF QEGSRELEASE LEAQLGQAEH RIRDLQSENQ RLKSEVDILK EKLEQQYAQS YKQISMLEDD LLQTRGIKEQ LHKYVRELEQ ANDDLERAKR ATITSLEDFE QRLNQAIIERN AFLESELDEK ELLVSVQRL KDEARDLRQE LAVRERTSDV TRMSAPSSPT LDIDKTDSAV QASLSLPATP VGKTMEHPFI GTKALTNGCG NGSPLTPSAR ISALNIVGDL LRKVGALSK LAACRNFAKD QAARKNYTTG NGNLINSNAT KFSHSLHTTY FDKTTMNGLD PGALSAAIAS PRAVSPPELL PLSV
Source	Mammalian cell
Target Names	ndel1b
Protein Names	Recommended name: Nuclear distribution protein nudE-like 1-B
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.