



Recombinant Mouse Sorting nexin-32 (Snx32)

Product Code	CSB-YP770506MO
Abbreviation	Snx32
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q80ZJ7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MEEQHQEAGN ESKPSSTSVD LQGDSPLQVE ISDAVSRDK VKFTVQTKSG LPHFAQSEFS VVRQHEEFIW LHDTYVENEE YAGLIIPPAP PRPDFEASRE KLQKLGEGNS SITREEFSKM KQELEAEYLA IFKKTVMAMHE VFLQRLAAHP TLRRDHNFSV FLEYSQDLSV REKNRKEVLG GLLRSIVRSA DEVLITGISG LKEVDDFFEHT ERTFLVEYHT RIRDTCQRAD RVMHSHKCLA DNYIPISAAL SSLGTQEVNQ LKRSFLKLAELFERLRKLEG RVASDEDLKL SDMLRYYMRD SQAAKDLLYR RLRALADYEN ANKALDKART RNREVRPAES RQQLCCQRFE RLSDSAKQEL MDFKSRRVSS FRKNLIELAE LELKHAKAST LLLQNTLVAL KGEP</p>
Source	Yeast
Target Names	Snx32
Protein Names	Recommended name: Sorting nexin-32 Alternative name(s): Sorting nexin-6B
Expression Region	1-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>