



Recombinant Human Properdin (CFP)

Product Code	CSB-BP005291HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P27918
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DPV LCFTQYEESS GKCKGLLGGG VSVEDCCLNT AFAYQKRSSG LCQPCRSRW SLWSTWAPCS VTCSEGSQLR YRRCVGNWNGQ CSGKVAPGTL EWQLQACEDQ QCCPEMGGWS GWGPWEPCSV TCSKGTRTRR RACNHPAPKC GGHCPGQAQE SEACDTQQVC PTHGAWATWG PWTPCSASCH GGPHEPKETR SRKCSAPEPS QKPPGKPCPG LAYEQRRCTG LPCCPVAGGW GPWGPVSPCP VTCGLGQTME QRTCNHPVPQ HGGPFCAGDA TRTHICNTAV PCPVDGEWDS WGEWSPCIRR NMKSISCQEI PGQQSRGRTC RGRKFDGHRG AGQQQDIRHC YSIQHCPLKG SWSEWSTWGL CMPPCGPNPT RARQLCTPL LPKYPPTVSM VEGQGEKNVT FWGRPLPRCE ELQGQKLVVE EKRPLHVPA CKDPPEEEL
Source	Baculovirus
Target Names	CFP
Protein Names	Recommended name: Properdin Alternative name(s): Complement factor P
Expression Region	28-469
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.