



Recombinant Fructose-like phosphotransferase enzyme IIB component 1 (fryB)

Product Code	CSB-MP770719SZB
Abbreviation	fryB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q821A9
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MSKKLIALCA CPMGLAHTFM AAQVLEEA AV EAGYEVKIET QGADGIQNRL TAQDIAEATI IHSVAVTPE DNERFESRDV YEITLQDAIK NAAGIIKEIE EMIASEQQ
Source	Mammalian cell
Target Names	fryB
Protein Names	Recommended name: Fructose-like phosphotransferase enzyme IIB component 1 EC= 2.7.1.69 Alternative name(s): PTS system fructose-like EIIB component 1
Expression Region	1-108
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.