



Recombinant Fumarate hydratase class II (fumC)

Product Code	CSB-MP771067SZB
Abbreviation	fumC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q83ML8
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MNTVRSEKDS MGAIQVPADK LWGAQTQRSL EHFRISTEKM PTSLIHALAL TKRAAAKVNE DLDLLSEEKA SAIRQAADEV LAGQHDDEFP LAIWQTGSGT QSNMNMNEVL ANRASELLGG VRGMRKVHP NDDVNKSQSS NDVFPTAMHV AALLALRKQL IPQLKTLTQT LSEKSRAFAD IVKIGRTHLQ DATPLTLGQE ISGWVAMLEH NLKHIEYSLP HVAELALGGT AVGTGLNTHP EYARRVADEL AVITCAPFVT APNKFEALAT CDALVQAHA GA LKGLAASLMK IANDVRWLAS GPRCGIGEIS IPENEPGSSI MPGKVNPTQC EALTMLCCQV MGNDVAINMG GASGNFELNV FRPMVIHNFL QSVRLLADGM ESFNKHCavg IEPNRERINQ LLNESLMLVT ALNTHIGYDK AAEIAKKAHK EGLTLKAAAL ALGYLSEAEF DSWVRPEQMV GSMKAGR
Source	Mammalian cell
Target Names	fumC
Protein Names	Recommended name: Fumarate hydratase class II Short name= Fumarase C EC= 4.2.1.2
Expression Region	1-467
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.