



Recombinant Protein FdhE homolog (fdhE)

Product Code	CSB-YP771081SZB
Abbreviation	fdhE
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q83PE8
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	MSIRIIPQDE LGSSEKRTAD MIPPLLPRL KNLYNRRRAER LRELAENNPL GDYLRFAALI AHAQKVLYD HPLEMDLTAR IKEASAQGKP PLDIHVLPRD KHWQKLLMAL IAELKPEMSG PALAVIENLE KASTQELEDM ASALFASDFS SVSSDKAPFI WAALSLYWAQ MANLIPGKAR AEYGEQRQYC PVCGSMPVSS MVQIGTTQGL RYLHCNLCET EWHVVRVKCS NCEQSGKLHY WSLDDEQAAI KAESCDDCGT YLKILYQEKE PKVEAVADDL ASLVLDARME QEGYARSSIN PFLFPGEGE
Source	Yeast
Target Names	fdhE
Protein Names	Recommended name: Protein FdhE homolog
Expression Region	1-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>