



Recombinant Arabidopsis thaliana Probable E3 ubiquitin-protein ligase ARI15 (ARI15)

Product Code	CSB-BP771158DOA
Abbreviation	ARI15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84RQ8
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MEKLMTESGL KPVVIDSNQD LSRVYCGICS NIGDDDDYDGD AVVVDGDLIS TPFCSHKFCK ACWSKYLKKN FFSVEKNHTA ISCPDRDCRA AVGPETVEKL TVRDQAMYEL YILKSYREKY LGWKLKLCPA RGCNYVIEFH LASEDEEHSL NIVCLCGHIF CWRCMLES HR PVTENNASDW LSRDLEK LIE EVDKPSTVSW IDANTKPCPH CFIPVEIDGE RPWAQFLTCV CSGRFCWKCF RSPETHGTSG SCLAPARSSN VGFNHNWRAK PGISCLDLWN ASQVNLVNAK YELEAFEESI IKKPSDLKEQ DVKVLREGLM LIVQCRQFLK WSCAYDYIHT EYDMAKREYL RFLQQNASGI VHSFSQSIKE ETEAKELTCG KLLSETT NIG NFFYHFIKTL REGLPEVQAE SYDNYGGPYW LCDRCTYGNS WFQRACKMCC DPTASKMDEL SD
Source	Baculovirus
Target Names	ARI15
Protein Names	Recommended name: Probable E3 ubiquitin-protein ligase ARI15 EC= 6.3.2.- Alternative name(s): ARIADNE-like protein ARI15 Protein ariadne homolog 15
Expression Region	1-452
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.